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Instructions for use

DIARFLEX® ADVANCED stool management system

EN

General information

Please read the instructions for use carefully each time before you use the products. The products must only be used by doctors, individuals instructed in their proper use and individuals entrusted with patient care.

Non-sterile: The DIARFLEX® ADVANCED stool management system is delivered non-sterile and is predominantly made from medical-grade silicone. All of the components of the system are free of latex.

1. Intended use

The DIARFLEX® ADVANCED stool management system is intended for the drainage of liquid and semi-liquid stool and ensures its secure collection in bedridden or immobile and incontinent patients (incl. irrigation and endoluminal administration of medication on a doctor's orders as well as sampling).

2. Product description

The DIARFLEX® ADVANCED stool management system is used for examining potentially infectious stool, contributes to reducing odours and preserving patient dignity and boosts the efficiency of care.

The DIARFLEX® ADVANCED stool management system comprises a silicone catheter tube with an odour-neutral coating, a collection bag with an activated charcoal filter and 2 replacement bags (each with a volume of 1.5 litres), a syringe, two hook-and-loop fasteners for affixing the bags and an additional bag holder.

A blue connector marked "45 ml" is used to fill/empty the low pressure balloon and a second connector marked "Irrig." provides access to the irrigation lumen. The enclosed syringe with the Luer lock connector can be used both for filling/emptying the balloon and for irrigating the system.

The central lumen of the silicone catheter tube includes a sampling port, which can be used to irrigate the catheter and take stool samples.

3. Period of use

Intended for a period of use of up to 29 consecutive days.

4. Indications

- Bedridden, immobile and incontinent patients
- Requirement for endoluminal medicament administration on a doctor's orders

5. Contraindications

- Do not use in patients with known sensitivity or allergies to materials used in this product.
- Do not use if the patient's distal rectum cannot accommodate the inflated volume of the low pressure balloon, or if the distal rectum/anal canal is severely strictured (anal stenosis).
- Do not use in patients with firm stools.
- Do not use in patients who have recently (under 6 weeks) undergone a rectal anastomosis or an anal or sphincter reconstruction.
- Do not use in patients with compromised rectal wall integrity (e.g. ischaemic proctitis, mucous membrane ulceration).
- Not suitable for paediatric use.

- Do not use in patients with pronounced haemorrhoids and corresponding symptoms.
- Do not use in patients with a suspected or diagnosed rectal/anal tumour.

6. Side effects

The following adverse events may be associated with the use of the stool management system:

- Perforation
- Pressure necrosis
- Loss of sphincter tone
- Obstruction
- Infection
- Excessive leakage of faecal matter

7. Warnings /Safety precautions

1. The enclosed warning labels for the filling lumen and irrigation lumen must be applied to the corresponding tubes before the system is used (blue warning label on blue filling tube, white label on white irrigation tube).
2. Only fill the low pressure balloon with the quantity of water specified in the Application section (45 ml, lukewarm). Do not use unpermitted quantities (overfilling) or types of liquid or hot liquids.
3. Rectal blood should be examined, in order to guarantee that there are no indications of pressure necrosis present. If pressure necrosis is identified, discontinue use and remove the catheter.
4. Do not connect mechanical pumping devices for the irrigation.
5. Do not irrigate patients with compromised intestinal wall integrity.
6. Excessive, prolonged traction on the catheter can lead to an erosion of the anal canal or the low pressure balloon can migrate into the anal canal.
7. Patients with too weak or too strong sphincter function may expel the catheter in normal use.
8. Extreme caution should be exercised in patients whose rectum has been modified by a radiation-induced stricture or whose rectum has been affected by radiation proctitis.
9. The patient must not lie on the drain tube or the connectors of the auxiliary lumens.
10. Do not sterilise and do not reuse.
11. DIARFLEX® ADVANCED is not indicated for use beyond 29 consecutive days.
12. Patients with severe tenesmus and patients who experience tenesmus or severe pain following insertion of the catheter may not tolerate the catheter. The catheter should be removed in such cases.

8. Application

8.1 Preparation


To achieve optimal performance of the catheter system, all formed stool should be evacuated from the colon prior to insertion of the system. Additionally, the rectum should be cleared of all stool prior to insertion of the catheter in all cases.

The preferred patient position for catheter insertion is the left lateral knee-chest position (the patient's clinical situation may dictate the use of an alternate position). The goal of patient positioning is to maximise sphincter relaxation in order to ease catheter insertion.

8.2. Insertion of the catheter


Screw the catheter tube onto the catheter bag.


Apply water-soluble lubricant generously to the deflated low pressure balloon and rectal opening.


Extract all the air from the low pressure balloon via the blue connector and then disconnect the syringe.


Grasp the lubricated catheter by the rear end of the low pressure balloon and introduce the catheter into the distal rectum using a rotational motion. Pay attention to the black reference ring indicating the maximum insertion depth

With a Luer Syringe, fill the the low pressure balloon with 40-45 ml of lukewarm water or saline solution via the blue connector and then disconnect the syringe.

Gently tug and release the low pressure balloon to position it and confirm that the catheter fits tension-free.

A catheter syringe (not included in the scope of supply) can be used to take a stool sample via the sampling port.

1. Moisten the tip of the catheter syringe.
2. Open the sampling port.
3. Remove an appropriately sized stool sample.
4. Close the sampling port again.

8.3. Care and monitoring of the system

- Check at regular intervals that the catheter and collection bag are positioned in such a way that the drain tube is not twisted, kinked or subject to external compression.
- Check that the catheter does not migrate into the colon (the black reference ring on the tube must always remain outside the anus).
- Check at regular intervals that waste is not accumulating in the catheter drain tube. (To resolve this issue, you can press the stool manually through the drain tube in the direction of the collection bag from outside with a guided movement.)


- Make sure and check at regular intervals the patient is not lying on the catheter drain tube or connectors of the auxiliary lumens.
- Flush the inside of the catheter drain tube at least twice daily (more frequently if required) and remove any accumulated faecal matter: To do so, inject water through the irrigation connector (marked with "Irrig.") a Luer syringe and then disconnect the syringe. Let the irrigation water flow into the collection bag.
- Replace the collection bag with a replacement bag when full.

8.4. Administration of medication

It is possible to administer medication rectally via the irrigation connector.


⚠ This may only be done on the attending doctor's orders and the doctor must check the dosage of medication.

- Perform a visual inspection to determine whether the irrigation line is clear and flush the line via the irrigation connector if necessary.
- Prepare Luer syringe filled with medication as per the doctor's orders.
- Position the enclosed catheter clamp at the black reference ring without closing it.


Insert the Luer syringe via the irrigation connector (connector marked „Irrig.”) and administer medication.

Flush the catheter immediately afterwards with tap water (50 ml, body temperature) via the irrigation connection (connector marked „Irrig.”) in order to transport the medication to the rectum.

Close the catheter clamp.


Once the medication exposure time specified by the doctor has elapsed, open and remove the catheter clamp.

8.5. Removal of the catheter

Use the Luer syringe to remove all the liquid from the low pressure balloon via the blue connector, then disconnect the syringe.

Apply water-soluble lubricant to the rectal opening, grasp the catheter just behind the entrance to the rectum and pull the low pressure balloon out carefully, applying even pressure.

⚠ Due to the elastic nature of the device, the removal procedure may result in the abrupt exiting of the catheter from the anus. Care should be taken to protect the caregiver from splatter.

9. Cleaning/ Disinfection

Do not clean or disinfect.

10. Storage/ Warehousing

The DIARFLEX® ADVANCED stool management system should be stored in a clean and dry environment protected from direct sunlight and heat.

11. Disposal

The devices must be disposed of in accordance with the respectively applicable national regulations governing waste materials.

12. Reporting

All serious events and incidents arising in connection with this medical device must be reported to the manufacturer and the responsible authorities.

13. Guaranteee

We guarantee that the delivered object is free from defects in accordance with the scope of our general terms and conditions (GTCs). You can find a copy of our GTCs among our business documentation and on our website at www.primed-halberstadt.de, etc.

Failure to comply with these terms will void the warranty immediately!

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Primed® Halberstadt Medizintechnik GmbH