ENG

GER

CareVis®

PROMECON

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
Manufactured exclusively for Promecon GmbH by:  
S&S Med CO., LTD.  
67 Bongseong-ro, Gunpo-si, Gyeonggi-do, 15848, Rep. of South Korea  
Tel +82 (31) 4723172. Fax +82 (31) 4723173

Easy Medical Device CommV  
Guldenberg 1 bus 2-4,  
2000 Antwerpen, Belgium  
Tel: +44 161 543 6588  
eurep@easymedicaldevice.com


Models / Modelle

## Models cont. / Modelle Forts.


Models with  
neuraxial patient  
connector in  
accordance with  
ISO 80369-6

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NaX models have a neuraxial patient connector as prescribed in the ISO 80369-6 for epidural, plexus block and wound infiltration anaesthesia procedures. The line (17) and the end cap (23) are yellow coloured to differentiate from the Luer-lock models.

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
The flow restrictor (4) is placed in-line.

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The total nominal flow rate is given by the flow-rate code ( **FF** ) and is written on the pump container. Each of the two distal infusion lines administers half of the total rate. **CareVis® DuO** flow rate labels on the restrictor of each individual infusion line show the nominal flow rate of this line only.

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## Priming the Infusion Line

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Filling time until reaching the patient connector (for pre-fixed flow-rate pumps only)

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OncO, PCA and TwiN pumps require longer than the specified time

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Email: info@promecon-medical.de, URL: www.promecon-medical.com

## Note

Please read the following information and instructions for use carefully. Failure to observe the instructions for use can endanger patient safety.

All serious incidents associated with the product should be reported to the manufacturer and as well to the competent authority of the EU Member State if the user and / or the patient is established within the EU.

## Models and Description

CareVis® is a single-use, portable infusion pump for continuous, on demand or combined continuous and on demand drug delivery that works independently from external energy sources like batteries or electricity. The required pump pressure to administer the solution is generated by an elastomeric balloon inside the infusion pump which is inflated with the solution. The flow rate, defined as the expelled fluid in ml/hr over the infusion period, is controlled by one or more flow restrictors incorporated into the infusion line.

CareVis® infusion pumps are available in different models as described below and, within these models, various volumes, flow rates, and other options can be selected.

The **CareVis®** models have a pre-fixed continuous flow-rate defined by the selected flow restrictor.

**CareVis® RuN** or **CareVis® DFO** are variants of the CareVis® models. The **RuN** models are quick-runners (voiding in 1/2 – 1hr) for the use with short-time antibiotic therapy. **DFO** models are designed for the infusion of deferoxamine, an iron chelating drug with a higher viscosity and a flow restrictor correspondingly adjusted to this viscosity (see under “Filling instructions” item 17).

The **CareVis® VariO** multi-rate models are also continuous pumps with the possibility to select seven different flow rates and to stop the flow with a detachable key (19).

The **CareVis® OncO** is a continuous pump with a wider but longer restrictor capillary (4a) suitable for use with 5-fluorouracil, which tends to precipitate easily.

The **CareVis® PCA** models offer the possibility to apply beside the continuous basal flow additional bolus doses via a patient-controlled analgesia (PCA) unit.

This pump is also provided without a continuous basal flow (**CareVis® BoluS** models).

The **CareVis® VariO Plus** models combine PCA function as described above with the multi-rate function of a CareVis® VariO within one unit (20). The two functions can be operated independently of each other.

The **CareVis® DuO** and **CareVis® TwiN** models have an infusion line branching into two infusion lines (21) as they may be used in wound infiltration anaesthesia with two catheters. The **CareVis® DuO** models have two flow restrictors, one in each line end, while the **CareVis® TwiN** models have one flow restrictor before the branching of the infusion line.

The **CareVis® NaX** models fulfil the criteria of the ISO 80369-6 for the neuraxial (epidural) use and for the use in plexus nerve block anaesthesia and in surgical wound (intra-lesional) anaesthesia. Instead of a Luer-lock patient connector they have the ISO 80369-6 patient connector (22) as designated for these indications. This is incompatible with the Luer patient catheters. Parts of the infusion line and the end cap (17, 23) are yellow coloured to signalise the ISO 80369-6 properties.

CareVis® single-use infusion pump system consists of an elastomeric balloon (1) inside a rigid, transparent container (2) and an infusion line (3) with a flow restrictor.

The protective transparent housing of the pump (2) provides UV blocking capabilities up to 390nm blocking UVB, UVC and most of UVA light.

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The flow restrictor is a capillary element (4) inside the patient connector (5) or in in-line position within the infusion line (4a) for the OncO models, or within the PCA (6), multi-rate units (18) or the combined PCA and multi-rate unit (20) respectively. The PCA unit (6, 20) has a bolus button to squeeze out an extra dose (7). When this bolus button of the PCA units (6, 20) is pressed, medication in a small bag beneath the bolus button is squeezed out and transferred to the patient. The re-filling time of this bag, generally called “lock-out time”, is pre-fixed. A flow restrictor within the line which re-fills the bag limits the flow into the bag at a fixed rate. More information about lock-out time is given under “Information for Patients” item 4.

The multi-rate units (18, 20) have a detachable key (19) to adjust the flow rate. A total of seven different flow rates can be set and the zero to stop the flow completely.

All CareVis® models are equipped with an in-line particle filter (8) of 1.2 micron, which is combined with a 0.2 micron pore size air filter (9). This last one forces trapped air out of the closed system. Its pore size of 0.2 micron is a barrier against bacteria from the outside.

The PCA pump models CareVis® PCA / BoluS / VariO and VariO Plus have the same filter system (8, 9) twice, one before the function units (6, 18, 20) and one after them.

All patient connectors are closed with a venting cap (14; 23 for NaX models), featured with a 0.2 µm bacteria filter. During filling and priming of pump and infusion line, the patient connector can be kept closed, as air but not fluid is escaping automatically through the venting cap (14, 23).

The flow rate accuracy of CareVis® infusion pumps fulfil the standards of ISO 28620:2020 for “Non-electrically driven portable infusion devices”: Infusion pumps with pre-fixed, continuous flow-rate: tolerance of ±15 % compared to the nominal flow rate and infusion pumps with adjustable flow rate: tolerance of ± 20 % compared to the nominal flow rate. In general, CareVis® infusion pump accuracy is within 10% and for adjustable flow rate within 15% resp.

## Indications of use

CareVis® single-use infusion pump is intended for the infusion of chemotherapy, antibiotic therapy, analgesic, or hydrating therapy and other general intravenous, intra-arterial, subcutaneous, epidural, perineural, and intra-lesional (in operative wounds) infusions.

CareVis® PCA, BoluS, and VariO or VariO Plus have been mainly designed for systemic or for regional pain management infusion use, and the corresponding NaX models for the use in plexus nerve blocks, the surgical wound infiltration anaesthesia, or in neuraxial (epidural) anaesthesia. Pain management infusion includes preoperative and postoperative surgery or chronic pain. The CareVis® TwiN and DuO models are for the use in surgical wound infiltration anaesthesia with two catheters.

CareVis® OncO is designed for the infusion chemotherapy with 5-fluorouracil. For other oncological infusions the regular CareVis® model is appropriate as well.

CareVis® RuN is developed for short-term antibiotic infusions and CareVis® DFO exclusively for the infusion of deferoxamine in iron eliminating treatment schemes.

All CareVis® models are suitable and safe for the use with children. However, the pumps should only be used for children and adult patients, who have fully comprehended and understood the function of the pump and, if selected, that of their optional PCA or multi-rate units, after having received adequate instructions by the healthcare provider.

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## Epidural Administration

Epidural administration of analgesics is limited to the use of indwelling catheters that are indicated specifically for epidural short-term or long-term administration of analgesics. In order to prevent the infusion of medicines not suitable for epidural use, no administration kits with side injection ports may be used. See also under "Precautions". It is strongly recommended that devices used for administration of medication via epidural routes be clearly differentiated from all other infusion devices. In conformity with the guidelines of the ISO 80369-6 the use of the special CareVis® NaX models is strongly recommended.

## Contraindications

Solutions with high viscosity or solid particles will cause slow flow or occlusion of the system and should not be used with CareVis® infusion pumps. As a result, also blood, blood products, proteins, and fluids containing lipids cannot be used.

**There is no alarm or alert when flow interruption occurs. Therefore, life-supporting medications, whose usage may cause serious injury or death due to stoppage or under-delivery, are contraindicated for infusion with CareVis® infusion pumps. The healthcare provider must ensure that a suitable alternative device for the infusion of critical medicines is used if an interruption of the infusion or changes in the infusion rate could result into serious health problems or even death.**

## Warnings

1.  CareVis® infusion pump is a single use product: DO NOT REUSE and DO NOT re-sterilise. Reuse or re-sterilisation may result in unpredictable flow of infusion and may result in infection of the patient.
2. Examine the patient before use to ensure that this product is suitable for the patient. The healthcare provider is responsible to ensure that the medical device is prepared and administered in accordance with the drug manufacturer's package leaflet. Physician is responsible for prescribing drug based on each patient's clinical status (such as age, body weight, disease state of patient, concomitant medications, etc.). Incorrect drug dosing or misuse may cause patient harm including drug toxicity.
3. It is the responsibility of the healthcare provider to ensure the patient is educated on the proper use of the medical device. It is the responsibility of the healthcare provider to modify Patient Guidelines provided with the pump as appropriate for your patients' clinical status and medication prescribed.
4. After filling CareVis® infusion pump, reconfirm the information and the flow rate settings prior to infusion, as incorrect flow rate or other incorrect parameters may fatally injure the patient.
5. **Epidural, perineural, or intralesional administrations of drugs other than those indicated for epidural, perineural, or intralesional use could result in serious irreversible injury and pain to the patient.**
6. Due to risk of ischemic injury, vasoconstrictors such as epinephrine are not recommended for infusions for the following routes of administration:
  - intraoperative site,
  - perineural, and
  - percutaneous.

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**7. There is no alarm or alert when flow interruption occurs. See under contraindications.**

8. If you use *CareVis® VariO* or *CareVis® VariO Plus* infusion pumps:

Flow rate is adjustable. To reduce potential adverse events, medication dosage should be calculated on maximum possible flow rate of the selected device.

The amount of medication over the therapeutic period and delivery time can vary by as much as 20%. Take this variance into consideration when determining medication delivery.

Regardless of the prescribed flow rate, only fill the pump with medication dosage that is appropriate to administer at the maximum possible flow rate of the device.

9. If you use *CareVis® PCA*, *CareVis® BoluS*, or *CareVis® VariO Plus* infusion pumps:

To reduce potential adverse events, medication dosage should be based on maximum applicable volume given by maximally using the bolus button.

**10. A parallel gravity infusion into the same patient catheter is strongly discouraged.** In case of a temporary catheter blockage, pumped drug may accumulate within the gravity line and suddenly released to the patient when the blockage opens.

### **Precautions**

CareVis® infusion pump will work correctly as long as the following precautions are respected.

1. Read the Instructions for Use carefully and don't use CareVis® infusion pump if you have not understood all details.

2. CareVis® infusion pump should be protected from sunlight and UV rays. Keep dry and store at temperatures between 4°C and 50°C.

3. Do not use CareVis® infusion pump if the packing is open or defective. Do not re-sterilise.

4. An aseptic technique must be used throughout the filling procedure.

5. The restrictor element (4, 4a) and the function units (6, 18, 20) may deteriorate or leak if placed in contact with disinfectants like alcohol.

6. Make sure the infusion line (3) is completely primed before use. Otherwise, air may be introduced into the patient.

7. The air-filter hole (9) on the air and particle filter (8) must not be covered and must be kept dry and clean. Solvents like alcohol may damage its seal against the infusion fluid. Equally, lipophilic (fat-soluble) drug solutions like those containing cremophor should not be used with the CareVis® infusion pumps.

8. Always supervise children and pets near the stored or in-use device.

9. Check before administration if the solution is clear. Do not use if opacities are visible.

10. During infusion, visually check the progress indicator (10) while the elastomeric balloon (1) gradually empties and ensure that infusion is duly proceeding.

11. The pump progress level indicator (10) on the protecting case does not represent the exact volume within the elastomeric balloon in ml. It is only intended to follow the pump progress by the patient or its caregiver. Infusion is complete when pump balloon is no longer inflated over the central piston.

12. No elements should be introduced through the air-inlet hole (11) of the container (2).

13. CareVis® infusion pump is allowed to come into contact with water, but it must be ensured that no appreciable amount of liquid enters the container (2) through the air-inlet hole (11).

14. Infusion should be interrupted by clamping (12) the line (3) and the medical caregiver informed if the pump body or its elements break or show leakage. Check the restrictor element (4, 4a), the function units (6, 18, 20), the air and particle filter (8) or any section of the infusion line.

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15. A residual volume of solution (around 1ml - 2ml depending on the pump size) will remain in the CareVis® infusion pump after the infusion is completed. See table “Pump Volumes and Volume Codes (V)” above.
16. At the end of infusion, dispose the device following the local standards on potentially toxic products.
17. When delivering drug solution through a catheter system, refer to the instructions for use provided by its manufacturer. The length, inner diameter, and position of catheter to be used can affect nominal flow rate. The use of catheters with small inner diameter (i.e. 22 or 23 Gauge or narrower) is strongly discouraged as it will slow down the flow.
18. It has been reported that there may be an association between the use of intra-articular catheters and/or the drugs (especially bupivacaine) infused intra-articularly and the development of chondrolysis. Avoid placing the catheter in joint spaces.
19. Avoid placing the catheter in the distal end of extremities (such as fingers, toes, nose, ears, penis, etc.) where fluid may build up as this may lead to ischemic injury or necrosis. Equally avoid tight wrappings which can limit blood supply or fluid diffusion.
20. **Note:** If pump did not perform as expected, **do not discard**. Contact your device provider or Promecon for product return instructions and helps so that the pump function can be investigated.

### Information on correct use

The accuracy of the flow rate depends on (1) filling volume, (2) temperature, (3) drug viscosity, (4) height difference between pump body and distal flow restrictor, and (5) used catheter size (see instructions for use).

1. CareVis® infusion pump is calibrated under the following conditions:
  - The pump filled at nominal volume as underfilling and overfilling, though permitted, will have influence on flow rate. See under “filling instructions” item 8.
  - Used infusion solution viscosity is that of normal saline (0.9% NaCl solution).
  - The temperature of the solution within the pump container (2) is at room temperature (23°C (73.5°F)).
  - The temperature of the flow restrictor (4, 4a) of CareVis®, CareVis® OncO, CareVis® RuN, CareVis® DFO and CareVis® DuO (including the NaX versions) is at 32°C (89.5°F), equivalent to the temperature when the restrictor is in contact with the skin.
  - The temperature of all other pumps (PCA, VariO, VariO Plus) including their function units (bolus or multi-rate) is at room temperature (23°C (73.5°F)).
  - The container (2) and the patient connector (5) are at the same level, with free output.
2. A variation of **viscosity** has an indirectly proportional effect on the nominal flow rate; an increase/decrease in viscosity slows down/increases the nominal flow rate.
3. A variation of **temperature** has a directly proportional effect on the actual flow rate. Flow rate will decrease or increase approximately 2% per 1°C (1.8°F) decrease or increase in temperature. Keeping the restrictor (4, 4a) of CareVis®, CareVis® OncO, CareVis® RuN, CareVis® DFO and CareVis® DuO (including the NaX versions) at direct contact with the skin of the patient helps to stabilise the operating temperature of 32°C (89.5°F). For all other pumps (PCA, VariO, VariO Plus) the patient or his caregiver has to care that the function units are held around room temperature and protected from sun light or cold weather.
4. Lifting the pump (2) above the distal connector (5) increases the **pressure gradient** and the actual flow rate, holding it below has the opposite effect.

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## Preparing the medication

1. The medication and its concentration to be infused must be prescribed by the patient's physician.
2. To prepare the solution, follow the specific instructions of the manufacturer for the medication in question, using an aseptic technique.
3. Choose the most suitable CareVis® infusion pump model for the volume and the flow rate as indicated on the identification label (13). This applies also to bolus volume and lock-out time if a PCA pump is used and to the flow-rate range if a multi-rate model is selected.

## Filling instructions (read carefully before filling)

1. An aseptic technique must be used throughout the procedure.
2. Do not remove the CareVis® infusion pump from its packaging until ready for use. Once the packaging has been opened, check that all the components are in operative condition and that the patient connector at the end of the infusion line is properly closed with the pre-installed venting cap (14, 23).
3. Remove the paper tape from the infusion line to unwind it.
4. Remove the sterility protection cap (15) from the filling port (16) and retain it aseptically for later use.
5. To determine the volume of solution to be infused, consider the residual volume (see table "Pump Volumes and Volume Codes (V)" above). Calculate the filling volume by multiplying the desired infusion time in hours by the nominal flow rate (ml/hr). Then, add the residual volume (see above under table "Pump Volumes and Volume Codes (V)"). An example: to prepare medication with a flow rate of 4ml/hr for 22 hours, the filling volume should be 89ml (4ml/hr x 22hrs + 1ml). The chosen infusion pump is an M40 (100ml nominal volume and 4 ml/hr flow rate). If a higher dilution of the drug and thus a larger infusion volume is desired, a larger pump with a higher flow rate is selected; in the same example an L100 pump filled with 222 ml and a flow rate of 10 ml/hr (22 hours x 10 ml/hr) + 2 ml residual volume).
6. Fill a Luer-lock syringe or a different filling device with the solution. Make sure that the syringe or filling device is completely drained of air. Attach the syringe or the line of a filling device to the filling port (16) by turning it gently and locking to the port. Do not over-tighten the syringe or the filling device in order not to damage the filling port. A needle or other sharp instrument may not be used. Fill the elastomeric balloon (1) with the calculated volume, if needed by repeating the operation. The solution should be introduced at a constant speed, without sudden spurts. Higher operating force may be required when starting the filling procedure of the elastomeric balloon (1). During filling of the elastomeric balloon, the medication flows throughout the infusion line (3) and removes automatically air from the infusion line through the in-line air filter (8) and through the filter of venting cap at the patient connector (14, 23). Both have hydrophobic (water-tight) filter membranes, which prevent leakage of the medication. The time required to prime the line should not exceed the time indicated in the table "Priming the Infusion Line".
7. The filling procedure may vary if **cytotoxic drugs** or drugs with tendency to precipitations like 5-fluorouracil are used. When cytotoxic drug is used, it is recommended to prime the patient line with normal saline first and fill up the active ingredient into the infusion pump immediately after. Initial priming of the infusion line with normal saline prevents tissue damage in case of initial opening extravasation. It also may impede precipitates within the infusion-line resulting in a decreasing drug flow rate. Therefore, small amounts of normal saline shall be filled first. When the solution has reached the patient connector (5) – this may take a while as indicated in the table "Priming the Infusion Line" – close the line with the clamp (12) to prevent osmotic immixture at its distal part. **Do not connect to the patient** if the distal parts of the infusion line (17) still contains air.

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8. The infusion pumps can be filled partly to its nominal volume but may cause an **increased flow-rate** than the nominal flow-rate, i.e. approx. 10% with the elastomeric balloon filled at 70% of nominal volume. Underfilling below 60% is not recommended.
- Overfilling the pump up to its indicated maximum filling volume may correspondingly **reduce the flow rate**.
9. The maximum volume of the CareVis® infusion pump as indicated in the “Pump Volume and Volumes Codes (V)” table above should not be exceeded.
10. Disconnect the syringe or the filling device from the filling port and close it (16) with the sterility protection cap (15) provided. Do not over-tighten the cap.
11. Fill out the patient label provided with all necessary information about the patient and the drug and affix onto the container (2).
12. Clamp (12) the infusion line near the patient connector if not done before.
13. Air bubbles may be trapped and visible inside the elastomeric balloon. This is normal and therefore do not try to remove trapped air bubble inside the balloon, as remaining air will be discharge by the air filter during the infusion.
14. **If you use CareVis® PCA**, the priming time is that of the continuous basal flow rate indicated in the “Priming the Infusion Line” table. In the meantime, the bolus bag fills within the designated filling time (lock-out time). The filling of the bag may expell small bubbles into the line distal of the bolus unit with the first pressing of the button; these bubbles will be eliminated by the second in-line air filter.
15. **If you use CareVis® BoluS** without a continuous basal flow (bolus-only models), the distal parts of the infusion line after the PCA button unit (17) **will not be primed** unless the bolus button is pressed for the first time after awaiting the selected lock-out time. The bolus bag beneath the bolus button will fill with some delay, which corresponds to the lock-out time, which is equivalent to the filling time of this bolus bag. Do not connect to the patient if the distal part of the infusion line (17) is not primed with the solution. Air will be pressed out completely by using the bolus button (7) having the slow filling time of the bolus bag in mind. Depending on the bolus volume selected, you may have to wait several times for the bolus reservoir to fill and then press again until the infusion line is completely filled. Have the lock-out time in mind before you repeat. When the line is completely filled, clamp (12) the line.
16. **If you use CareVis® VariO or CareVis® VariO Plus** it is recommended to select the highest possible flow rate (3, 5, 7 or 14 ml/hr depending on selected model) for efficient priming. After successfully priming the infusion line up to the patient connector, set the flow-rate to the flow rate chosen for the patient to start the infusion or to zero. In both cases clamp (12) the infusion line near the patient connector (17). **Do not fill VariO pumps with the flow rate set to zero** as the distal part of the line beyond the VariO unit will not fill. **Do not connect to the patient** if the distal parts of the infusion line (17) still contains air.
17. **If you use CareVis® DFO** adjustments of the concentration (by the filling volume) may be necessary depending on the total dose of deferoxamine as the restrictor is calibrated for a concentration of 2 g deferoxamine in 96 ml normal saline.

## Application to patient

1. Check the patient label for correct patient name and correct drug and drug concentration data. Check the parameters of the selected infusion pump model before use and check if these fulfil the indication and needs of the patient: The filling volume and infusion rate are labelled on the pump body, bolus information on the bolus unit and variable rate information on the multi-rate unit.

**Important: Verify that the route of administration is correct.**

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2. Before starting infusion, check that the filled device is at room temperature. If the unit has been refrigerated, it should be equilibrated at room temperature before use, which may take 8 to 18 hours depending on pump volume.

Filling Volume (ml)	60/100	180	275	450	600
Refrigerator to Room Temp (hrs)	8	10	12	15	18

Do not microwave or use an external heat source. Dew may form on the inner surface of the case (2) during this time but does not adversely affect the normal use of the unit. Dew may also form with minor changes of temperature.

3. Before starting infusion, remove the blue (14) (yellow (23) in NaX models) infusion line end cap from distal patient connector (5), open the clamp (12) and visually verify that the solution has reached the tip of distal patient connector. The device is vented when two drops of drug solution emerge at the patient connector.

4. Close the clamp again and connect the distal connector firmly to the device used for the infusion to patient, open the clamp and start infusion.

5. For CareVis<sup>®</sup> infusion pumps with pre-fixed basal flow (CareVis<sup>®</sup>, CareVis<sup>®</sup> OncO, CareVis<sup>®</sup> RuN, CareVis<sup>®</sup> DFO and CareVis<sup>®</sup> DuO (including the NaX versions)), attach the flow restrictor (4, 4a), to the patient's skin using a tape or equivalent dressing.

For CareVis<sup>®</sup> infusion pumps with a bolus, a multi-rate or a combined unit (CareVis<sup>®</sup> PCA, CareVis<sup>®</sup> BoluS, CareVis<sup>®</sup> VariO, CareVis<sup>®</sup> VariO Plus, and CareVis<sup>®</sup> VariO TwiN (including the NaX versions)), please make sure these function units are operated at room-temperature condition.

6. For all CareVis<sup>®</sup> infusion pumps, ensure that the air and particle filter (8) is kept dry and the air-filter hole (9) is not obstructed. Do not tape the filter.

7. **In case a CareVis<sup>®</sup> PCA or a CareVis<sup>®</sup> BoluS or a CareVis<sup>®</sup> VariO Plus is used**, instruct the patient carefully about its use and the clinical condition when to press the bolus button or better not.

8. **In case a CareVis<sup>®</sup> VariO or CareVis<sup>®</sup> VariO Plus is used**, set the initial flow rate with the setting key and detach it (19). If it is decided to leave the key and the selection of the subsequent flow rates to the patient, **instruct the patient carefully when and how to do it and about the risks of a too high flow rate of pain drug**. In general, the setting of the flow rates is done by the nursing personnel or the physician and the key is not left with the patient.

9. The infusion is finished when the level indicator (10) is close to 0 and the elastomeric balloon (1) is empty and fully deflated over the central piston.

10. Visually inspect the device once the infusion has finished. Clamp (12) the infusion line and disconnect the patient connector from the patient (catheter). To avoid possible handling risks, connect the end of the output connection (5) to the Luer-Lock filling port (16) (not possible with NaX models).

11. Dispose of the device following the centre's protocol or following the local norms.

## Information for patients

Instruct the patient about the following:

1. The level indicator on the infusion pump container (10) is used to periodically monitor the medication infusion process. If the indicator does not change throughout the infusion, the healthcare provider should be informed.
2. Do not bring soap or alcoholic solutions in contact with the filter (8), the restrictor element (4, 4a), the PCA unit (6), or the multi-rate unit (18, 20), as this may cause leakage from the air-eliminating vent or the joints.

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3. Do not disconnect the device before infusion is completed unless otherwise directed by the healthcare provider. The infusion is finished when the level indicator (10) is close to zero and the elastomeric balloon (1) is empty and fully deflated (tightly attached to the underlying piston).
4. **In case the *CareVis® PCA* or a *CareVis® BoluS* or a *CareVis® VariO Plus* is used**, the PCA bolus unit provides an extra dose of pain medication if the degree of pain requires it. The bolus bag (7) should always be discharged completely by pressing the bolus button all the way down and then been released. Follow the instructions of your doctor or the nursing personnel about the use and the frequency of use of the bolus button. The **lock-out time** of the *CareVis® PCA*, the *CareVis® BoluS*, or the *CareVis® VariO Plus* is the required time to fill the bolus bag, which is located beneath the bolus button (7). Pressing the bolus button before the designated lock-out time will cause a proportionally reduced bolus dose as the bolus bag would not be completely filled. Pressing the bolus button repeatedly within the lock-out time will not deliver more medication than the designated total bolus volume of medication selected by your doctor for you. If you believe that you need more pain medication than given within the lock-out time, contact your healthcare provider (doctor or nurse) early. The tolerance of maximal doses may vary from patient to patient.
5. **If you are attached to a *CareVis® VariO* or *CareVis® VariO Plus*** the flow rate can be changed with the help of the detached key (19). Usually this is done when the pain level changes (either way) either by the nursing personnel or the doctor or by the patient if the patient is instructed to do so and correspondingly informed.
6. In order to maintain pump accuracy, you should wear the pump on body height and when in bed put it down on the same height as your bed, for example on your bed or your bedside table. Do not hang the pump above you or put it on the floor. **Do not put the pump below your blanket.**
7. If you receive a pain medication with the *CareVis®* infusion pump notify your doctor when following signs or symptoms appear:
- Dizziness, drowsiness, straining vision, buzzing or ringing in the ears
  - Metal taste in the mouth
  - Tingling and / or numbness around the mouth, the finger ends or toes
  - Sleepiness or confusion
  - Breathing difficulties
- Your doctor may inform you about other reasons to contact him during infusion.
8. **If you use a *CareVis®*, *CareVis® OncO*, *CareVis® RuN*, *CareVis® DFO* or *CareVis® DuO*** pay attention that the end piece of the infusion line is attached to your skin in order to guarantee skin contact at all times.
9. **If a pump with a functional unit is used (*CareVis® PCA*, *CareVis® BoluS*, *CareVis® VariO* or *CareVis® VariO Plus*)**, this unit should be kept as close as possible to room temperature, i.e. not under clothing or the blanket.
10. The pump, the filter and the wound site should be protected from water. Do not drown the pump into the water. Single water splashes are harmless.

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