

MULTIFLOW[®]



The Multiflow System

Is specially designed for postoperative pain management. It offers a wide range of flow rates that can be adjusted to patient/physician needs, from 1ml/h to 7ml/h in discrete steps of 1ml/h, and from 2ml/h to 14ml/h in discrete steps of 2ml/h.



Leventon

A Werfen Company

MULTIFLOW[®] Advantages

1. Therapy Flexibility

Freedom to adjust the infusion rates allowing the physician to customise the patient's pain management protocol.

2. Controlled Access

Safe lock with removable key provides flow control only to selected individuals.

3. Responsive Pain Management

Amount of drug infusion can be adapted to changing patient's needs.

4. Effective Pain Management

Helps to reduce hospital stay and improve recovery time.

5. Robust Container

Rigid plastic housing with a clear fluid level Indicator.



100ml	250ml	400ml	600ml
<p>100 ml 1 - 7ml/h 100-124</p>	<p>250ml 1 - 7ml/h 250-124</p>		
<p>100 ml 2 - 14ml/h 100-248</p>	<p>250ml 2 - 14ml/h 250-248</p>	<p>400 ml 2 - 14ml/h 400-248</p>	<p>600 ml 2 - 14ml/h 600-248</p>

PL-259000346-1V07 06-2022



Distributed by:

Manufactured by:

ERMED AG

Tel 052/687 40 80 Fax 052/687 40 89
www.erved.ch info@erved.ch
 Zunderst Wyler 14 8226 Schleithelm

Leventon
 A Werfen Company

Leventon, S.A.U.
 C. Newton, 18-24
 08635 Sant Esteve Sesroviros (Barcelona). Spain
 Tel: +34 938 176 300
 Fax: +34 938 176 301
leventon.es