

Epidural Needle



tuohy

- Perfect angular shape for excellent application



Sizes	reference no.	product name	gauge	color	length mm	length inch
	TEP1690	Epidural Needle 16G - 90 mm	16G		90mm	3½"
	TEP1880	Epidural Needle 18G - 80 mm	18G		80mm	3"
	TEP1890	Epidural Needle 18G - 90 mm	18G		90mm	3½"
	TEP18100	Epidural Needle 18G - 100 mm	18G		100 mm	4"

EPIFIX[®]

Epidural Anesthesia Sets



Catheter Guide

Epidural Catheter



Epidural LOR Syringe



EPIjet Automatic LOR Syringe



Epidural Needle



Catheter Connector



Catheter Fixing Device



Epidural Filter

Filter Fixing Device



Epidural Anesthesia Sets

ORDER TABLE	REFERENCE NO	TUOHY NEEDLE				EPISAFE STANDARD 20G 3 LATERAL HOLES	EPISOFT SOFT TIP 20G 3 LATERAL HOLES	EPIUNO ATRAUMATIC 20G 6 LATERAL HOLES	EPIUNO SOFT ATRAUMATIC SOFT TIP 20G 6 LATERAL HOLES	EPILOL LOR SYRINGE 10 ML.	EPIJET AUTOMATIC LOR SYRINGE 10 ML.	EPIDURAL FILTER
		16 G-90mm	18 G-90mm	18 G-100mm	18 G-80mm							
EPIFIX SET 16G	TEPIS16903	*				*				*		*
	TEPIX16903	*				*					*	*
	TSEPIS16903	*					*			*		*
	TSEPIX16903	*					*				*	*
	TEPISU16906	*						*		*		*
	TEPIXU16906	*						*			*	*
	TSEPISU16906	*							*	*		*
	TSEPIXU16906	*							*		*	*
EPIFIX SET 18G	TEPIS18903		*			*				*		*
	TEPIX18803				*	*				*		*
	TEPIX18903		*			*					*	*
	TEPIX18803				*	*					*	*
	TSEPIS18100			*			*			*		*
	TSEPIS18903		*				*			*		*
	TSEPIS18803				*		*			*		*
	TSEPIX18903		*				*				*	*
	TSEPIX18803				*		*				*	*
	TEPISU18906		*					*		*		*
	TEPISU18806				*			*		*		*
	TEPIXU18906		*					*			*	*
	TEPIXU18806				*			*			*	*
	TSEPISU18906		*						*	*		*
	TSEPISU18806				*				*	*		*
	TSEPIXU18906		*						*		*	*
	TSEPIXU18806				*				*		*	*
EPIFIX MINI A SET	TEPIMA16903	*				*						*
	TSEPIMA16903	*					*					*
	TEPIMAU16906	*						*				*
	TSEPIMAU16906	*							*			*
	TEPIMA18903		*			*						*
	TSEPIMA18903		*				*					*
	TEPIMAU18906		*					*				*
	TSEPIMAU18906		*						*			*
EPIFIX MINI B SET	TEPIMB3					*						*
	TSEPIMB3						*					*
	TEPIMBU6							*				*
	TSEPIMBU6								*			*
EPIFIX MINI C SET	TEPIMC16903	*				*						
	TSEPIMC16903	*					*					
	TEPIMCU16906	*						*				
	TSEPIMCU16906	*							*			
	TEPIMC18903		*			*						
	TSEPIMC18903		*				*					
	TEPIMCU18906		*					*				
	TSEPIMCU18906		*						*			

Epidural Set Content

- Epidural Needle "Tuohy" (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- Epidural Filter
- Catheter Fixing Device
- Filter Fixing Device
- Catheter Connector

Epidural Mini Set Content

- Epidural Needle "Tuohy"(See the Table)
- Epidural Catheter (See the Table)
- Epidural Filter(See the Table)
- Catheter Connector
- Catheter Fixing Device

COMBIFIX[®]

Combined Spinal-Epidural Anesthesia Sets (CSE)

twist - lock Lock System



Epidural LOR Syringe



EPIject Automatic LOR Syringe



Epidural Catheter

Catheter Guide



Catheter Connector



Epidural Filter



Catheter Fixing Device



Filtered Aspiration Cannula
■ It enables filtration of the anesthetic fluid aspirated from ampule

COMBIFIX®

Combined Spinal-Epidural Anesthesia Sets (CSE)

ORDER TABLE	REFERENCE No.	TUOHY NEEDLE			PENCIL POINT SPINAL NEEDLE 27G	QUINCKE TIP SPINAL NEEDLE 27G	EPISAFE STANDARD 20G 3 LATERAL HOLES	EPISOFT SOFT TIP 20G 3 LATERAL HOLES	EPIUNO ATRAUMATIC 20G 6 LATERAL HOLES	EPIUNO SOFT ATRAUMATIC 20G 6 LATERAL HOLES	EPILOL LOR SYRINGE 10 ML	EPIJET AUTOMATIC LOR SYRINGE 10 ML
		16 G-90mm	18 G-90mm	18 G-90mm Back Eye								
COMBIFIX SET 16G	TSWS163	*			*		*				*	
	TSWX163	*			*		*					*
	TSWSS163	*			*			*			*	
	TSWXS163	*			*			*				*
	TSWSU166	*			*				*		*	
	TSWXU166	*			*				*			*
	TSWSUS166	*			*				*		*	
COMBIFIX SET 18G	TSWS183		*		*		*				*	
	TSWSS183		*		*			*			*	
	TSWXS183		*		*			*				*
	TSWSU186		*		*				*		*	
	TSWXU186		*		*				*			*
	TSWSUS186		*		*					*	*	
	TSWXUS186		*		*					*		*
	TSWS183Q		*			*	*				*	
	TSWX183Q		*			*	*					*
	TSWSS183Q		*			*		*			*	
	TSWXS183Q		*			*		*				*
	TSWSU186Q		*			*			*		*	
	TSWXU186Q		*			*			*			*
	TSWSUS186Q		*			*				*	*	
	TSWXUS186Q		*			*				*		*
COMBIFIX BACK EYE	TSWSS18D3			*	*			*			*	
	TSWXS18D3			*	*			*				*
	TSWSS18D3Q			*		*		*			*	
	TSWXS18D3Q			*		*		*				*
COMBIFIX 18G Mini SET	TSWM16	*			*						*	
	TSWXM16	*			*							*
	TSWM18		*		*						*	
	TSWXM18		*		*							*

Combined Set Content

- Epidural Needle "Tuohy" (See the Table)
- Spinal Needle (See the Table)
- Epidural Catheter (See the Table)
- 10 ml Epidural LOR Syringe (See the Table)
- Epidural Filter
- Filtered Aspiration Cannula
- Catheter Fixing Device
- Filter Fixing Device
- Catheter Connector

Combined Mini Set Content

- Epidural Needle "Tuohy" (See the Table)
- Spinal Needle (See the Table)
- Epidural LOR Syringe (See the Table)
- Filtered Aspiration Cannula