

Spinal Anesthesia Needles



*Optional 25G-26G-27G-29G

ATRAUCOM®

- Multi-angle grinded tip with special bending
- Extra thin wall
- Atraumatic

PenPen® pencil point


- Perfect conic elliptical bevel enables easy access to desired location, without cutting or causing harm


Filtered Aspiration Cannula
It provides filtration of the anesthetic fluid aspirated from the ampule.
*Optional


QPen® quincke tip


- Quincke bevel shaped and smoothened by using state of the art technology goes into tissues easily.
- While perfect match of stylet and needle bevel prevents tissues to go in the cannula, angular tip prevents tissue damage.
- Crystal clear needle hub enables monitoring fluid flow and helps to keep CSF loss at a minimum level.
- Thin needle wall enables rapid and easy flow.

Spinal Anesthesia Needles

Sizes	reference no.	product name	gauge	color	length mm	length inch	guide needle
QPen®  quincke	TSPQ1875	18 G – 75 mm Spinal Needle Quincke Tip	18G	pink	75mm	3"	
	TSPQ1890	18 G – 90 mm Spinal Needle Quincke Tip	18G	pink	90mm	3½"	
	TSPQ1990	19G – 90 mm Spinal Needle Quincke Tip	19G	orange	90mm	3½"	
	TSPQ2075	20 G – 75 mm Spinal Needle Quincke Tip	20G	yellow	75mm	3"	
	TSPQ2090	20 G – 90 mm Spinal Needle Quincke Tip	20G	yellow	90mm	3½"	
	TSPQ2190	21 G – 90 mm Spinal Needle Quincke Tip	21G	green	90mm	3½"	
	TSPQ2240	22 G – 40 mm Spinal Needle Quincke Tip	22G	black	40mm	1½"	
	TSPQ2275	22 G – 75 mm Spinal Needle Quincke Tip	22G	black	75mm	3"	
	TSPQ2290	22 G – 90 mm Spinal Needle Quincke Tip	22G	black	90mm	3½"	
	TSPQ22120	22 G – 120 mm Spinal Needle Quincke Tip	22G	black	120mm	4¾"	
	TSPQ2390	23G – 90 mm Spinal Needle Quincke Tip	23G	blue	90mm	3½"	
	TSPQ2490	24 G – 90 mm Spinal Needle Quincke Tip	24G	purple	90mm	3½"	
	TSPQ2575	25 G – 75 mm Spinal Needle Quincke Tip	25G	orange	75mm	3"	
	TSPQ2590	25 G – 90 mm Spinal Needle Quincke Tip	25G	orange	90mm	3½"	
	TSPQK2590	25 G – 90 mm Spinal Needle Quincke Tip, with Guide	25G	orange	90mm	3½"	20G
	TSPQK25103	25 G – 103 mm Spinal Needle Quincke Tip, with Guide	25G	orange	103mm	4"	20G
	TSPQK25120	25 G – 120 mm Spinal Needle Quincke Tip, with Guide	25G	orange	120mm	4¾"	20G
	TSPQ2690	26 G – 90 mm Spinal Needle Quincke Tip	26G	brown	90mm	3½"	
	TSPQK2690	26 G – 90 mm Spinal Needle Quincke Tip, with Guide	26G	brown	90mm	3½"	20G
	TSPQK26103	26G – 103 mm Spinal Needle Quincke Tip, with Guide	26G	brown	103mm	4"	20G
	TSPQK26120	26G – 120 mm Spinal Needle Quincke Tip, with Guide	26G	brown	120mm	4¾"	20G
	TSPQ2790	27 G – 90 mm Spinal Needle Quincke Tip	27G	grey	90mm	3½"	
	TSPQK2790	27 G – 90 mm Spinal Needle Quincke Tip, with Guide	27G	grey	90mm	3½"	22G
	TSPQK27103	27 G – 103 mm Spinal Needle Quincke Tip, with Guide	27G	grey	103mm	4"	22G
	TSPQK27120	27 G – 120 mm Spinal Needle Quincke Tip, with Guide	27G	grey	120mm	4¾"	22G
	TSPQK2990	29G – 90 mm Spinal Needle Quincke Tip, with Guide	29G	red	90mm	3½"	23G
	TSPQK29120	29G – 120 mm Spinal Needle Quincke Tip, with Guide	29G	red	120mm	4¾"	23G

Sizes	reference no.	product name	gauge	color	length mm	length inch	guide needle
PenPen®  pencil point	TSPK2290	22 G – 90 mm Spinal Needle Pencil Point	22G	black	90mm	3½"	
	TSPPK2290	22 G – 90 mm Spinal Needle Pencil Point, with Guide	22G	black	90mm	3½"	18G
	TSPK2390	23 G – 90 mm Spinal Needle Pencil Point	23G	blue	90mm	3½"	
	TSPK2490	24 G – 90 mm Spinal Needle Pencil Point	24G	purple	90mm	3½"	
	TSPK2590	25 G – 90 mm Spinal Needle Pencil Point	25G	orange	90mm	3½"	
	TSPPK2590	25 G – 90 mm Spinal Needle Pencil Point, with Guide	25G	orange	90mm	3½"	20G
	TSPPK25103	25 G – 103 mm Spinal Needle Pencil Point, with Guide	25G	orange	103mm	4"	20G
	TSPPK25120	25 G – 120 mm Spinal Needle Pencil Point, with Guide	25G	orange	120mm	4¾"	20G
	TSPK2690	26 G – 90 mm Spinal Needle Pencil Point	26G	brown	90mm	3½"	
	TSPPK2690	26 G – 90 mm Spinal Needle Pencil Point, with Guide	26G	brown	90mm	3½"	20G
	TSPPK26103	26 G – 103 mm Spinal Needle Pencil Point, with Guide	26G	brown	103mm	4"	20G
	TSPPK26120	26 G – 120 mm Spinal Needle Pencil Point, with Guide	26G	brown	120mm	4¾"	20G
	TSPK2790	27 G – 90 mm Spinal Needle Pencil Point	27G	grey	90mm	3½"	
	TSPPK2790	27 G – 90 mm Spinal Needle Pencil Point, with Guide	27G	grey	90mm	3½"	22G
	TSPPK27103	27 G – 103 mm Spinal Needle Pencil Point, with Guide	27G	grey	103mm	4"	22G
	TSPPK27120	27 G – 120 mm Spinal Needle Pencil Point, with Guide	27G	grey	120mm	4¾"	22G
	TSPK2990	29 G – 90 mm Spinal Needle Pencil Point, with Guide	29G	red	90mm	3½"	23G
	TSPPK2590	25 G – 90 mm Pencil Point, with Guide, Filtered Aspiration Cannula	25G	orange	90mm	3½"	20G
	TSPPK2690	26 G – 90 mm Pencil Point, with Guide, Filtered Aspiration Cannula	26G	brown	90mm	3½"	20G
	TSPPK2790	27 G – 90 mm Pencil Point, with Guide, Filtered Aspiration Cannula	27G	grey	90mm	3½"	22G

Sizes	reference no.	product name	gauge	color	length mm	length inch	guide needle
ATRAUCOM® 	TSPAK2590	25 G – 90 Spinal Needle, with Guide	25G	orange	90mm	3½"	20G
	TSPAK25100	25 G – 100 mm Spinal Needle, with Guide	25G	orange	100mm	4"	20G
	TSPA2650	26 G – 50 mm Spinal Needle	26G	brown	50mm	2"	
	TSPAK2688	26 G – 88 mm Spinal Needle, with Guide	26G	brown	88mm	3½"	20G
	TSPAK26100	26 G – 100 mm Spinal Needle, with Guide	26G	brown	100mm	4"	20G


Sizes	reference no.	product name	gauge	color	length mm	length inch
Guide Needle 	KIH18	18 G – 34 mm Guide Needle	18G	pink	34mm	1½"
	KIH20	20 G – 34 mm Guide Needle	20G	yellow	34mm	1½"
	KIH21	21 G – 34 mm Guide Needle	21G	green	34mm	1½"
	KIH22	22 G – 34 mm Guide Needle	22G	black	34mm	1½"
	KIH23	23 G – 34 mm Guide Needle	23G	blue	34mm	1½"

SPINOSET[®]

Spinal Anesthesia Sets




5 ml Syringe




3 ml Syringe



Medical Hydrophilic Swab



Spinal Needle



Guide Needle



Square Plaster



Waste Needle Box



Hypodermic Needle
(Local Anesthesia)



Sponge



Sterile Surgical Towel




















































Transparent Adhesive Surgical Drape



Filtered Aspiration Cannula
It enables the filtration of the
anesthetic fluid aspirated from
flacon or ampule

SPINOSET®

Spinal Anesthesia Sets

ORDER TABLE	REFERENCE NO.	GAUGE	COLOR				
				QUINCKE SPINAL NEEDLE	PENCIL POINT SPINAL NEEDLE	ATRAUCOM SPINAL NEEDLE	GUIDE NEEDLE
QUINCKE TIP SET	PSPQ1875	18G		*			
	PSPQ1890	18G		*			
	PSPQ1990	19G		*			
	PSPQ2075	20G		*			
	PSPQ2090	20G		*			
	PSPQ2190	21G		*			
	PSPQ2240	22G		*			
	PSPQ2275	22G		*			
	PSPQ2290	22G		*			
	PSPQ22120	22G		*			
	PSPQ2390	23G		*			
	PSPQ2490	24G		*			
	PSPQ2575	25G		*			
	PSPQK2590	25G		*			*
	PSPQK25103	25G		*			*
	PSPQK25120	25G		*			*
	PSPQK2690	26G		*			*
	PSPQK26103	26G		*			*
	PSPQK26120	26G		*			*
	PSPQK2790	27G		*			*
	PSPQK27103	27G		*			*
	PSPQK27120	27G		*			*
	PSPQK2990	29G		*			*
	PSPQK29120	29G		*			*
PENCIL POINT SET	PSPP2290	22G			*		
	PSPPK2290	22G			*		*
	PSPP2390	23G			*		
	PSPPK2390	23G			*		*
	PSPP2490	24G			*		
	PSPPK2490	24G			*		*
	PSPPK2590	25G			*		*
	PSPPK25103	25G			*		*
	PSPPK25120	25G			*		*
	PSPPK2690	26G			*		*
	PSPPK26103	26G			*		*
	PSPPK26120	26G			*		*
	PSPPK2790	27G			*		*
	PSPPK27103	27G			*		*
	PSPPK27120	27G			*		*
	PSPPK2990	29G			*		*
ATRAUCOM SET	PSPAK2590	25G				*	*
	PSPAK25100	25G				*	*
	PSPA2650	26G				*	
	PSPAK2688	26G				*	*
	PSPAK26100	26G				*	*