

EYE SAFETY TAPE

ANAESTHESIA - OPERATION ROOM



ERMED AG

Tel 052/687 40 80
www.armed.ch
Zunderst Wyler 14


Fax 052/687 40 89
info@armed.ch
8226 Schleithem

 **Pergo®**
Passion for Better

Intended Use

Eye Safety Tape is used for covering the eyes during surgeries while patient is under anaesthesia. During anaesthesia, most of the time the eyes are left open which result in dry eyes. Dry eyes can lead to eye injuries. Currently various gels are used to avoid drying on the eyes and nurses have to check if eyes dried out on constant basis. Current common application to avoid drying eyes is medical tape. This application harms skin, eyebrows and eyelashes. Medical tape is not transparent to see the position of the eyes and it does not fit the eye anatomy properly.


According to various studies, 59% of the eyes left open during operations and it is recommended to close and protect the eyes while patient is under anaesthesia.



05TD0012 Eye Safety Tape Adult

Key Features

- Coloured, ergonomic and non sticky handle.
- Designed to fit eye anatomy perfectly.
- Transparent material.
- Hypoallergic and medical grade material.
- Latex free.
- Patent protected.



05TD0020 Eye Safety Tape Pediatric

Advantages for Patients

- * Does not hurt the skin, eyebrows and eyelashes while removing (post operation).
- * Protects eyes from external splashes during the operation.
- * Eyes do not dry out even during long operations.
- * Medical grade material avoids any skin irritation (even on sensitive skin).
- * Reduces corneal abrasion risk.

Advantages for Healthcare Professionals

- * Transparent structure makes for easy pupil dilation control visible at any time without removing the Eye Safety Tape.
- * Coloured and non sticky handle allows easy and fast application and removal. Reduces operational risk and saves on application time (cost).